

## ROLLBOARD / PATIENT TRANSFER SLIDE ASM R200

\* Easy and Secure Patient Transfer

\* Lightweight and Practical To Use

\* Color Options

\* Carrierbag and Wall Holder



### Technical Features

- \* It has a flexible eva sheet inside it.
- \* The eva sheet is covered with polyester fabric.
- \* There is a slippery fabric that makes it easy to turn the plate around.
- \* The rollboard is easy to double up, it can be carried from the holding handle.
- \*Color Options: Red, Blue
- \*Sizes: When it's closed 47x85cm  
When it's open 47x172cm
- \*The Weight: 2.4 kg
- \*Board Thickness: 1.9cm
- \*Transfer Capacity: 160kg
- \*X-Ray Permeability
- \*The rollboard is easy to clean up
- \*Antibacterial

### Standart Equipment

- \* Rollboard
- \* Carrierbag
- \* Wall Holder

### Product image



### Wall Holder



### Carrierbag

