



# asm

ASM TIBBİ CİHAZLAR ve MÜH.  
HİZ. İNŞ. SAN. TİC. LTD. ŞTİ.

Karlıktepe Mah. Ayabakan Sok. No:5-7/B Kartal/İSTANBUL  
Tel.: +90 216 488 50 00 Faks: +90 0216 488 18 28  
info@asmmed.com www.asmmed.com

## ROLLBOARD / PATIENT TRANSFER SLIDE ASM R100

In health institutions, operating rooms, clinics, in short; the **Rollboard** is being used everywhere to transfer the patients to different locations when needed in a practical way . This product is manufactured using extremely durable materials, it is specially designed to ensure that the patients are being transferred with minimum jolt.

- \* The Rollboard can be perfectly used to transfer patients after their surgery because it causes minimum jolt.
- \* There is a special fabric on it. The patient can be easily put on the product and easily transferred to a different patient bed because of this fabric.
- \* The patient can be carried to the stretcher or the patient bed without being moved.
- \* The Rollboard is so designed that the patient can be transferred with just a single person.
- \* “Custom manufacturing in desired color” For example: The product in the operating room and services can be used in different colors so the slides can be prevented from being switched



WALL HOLDER



- \* 150 kg carrying capacity
- \* it is 1800 mm long, 500 mm wide,  
The height is only 20 mm.
- \* Antibacterial and easy to be cleaned
- \* It's made of x-ray permeable material

- \* Wall Holder
- \* It's easy to double up
- \* It's easy to carry with the shoulder straps
- \* Strong and changeable cover
- \* Easy and safe transfer