

Electronic personal scale

M320600-01

ADE



The perfect choice for use in all medical fields. Heavy duty electronic column scale with an extra low level slip resistant platform. The transport castors assure an easy movability. Straight forward BMI, Auto-Hold and mother-child operation gives fast and clear results of the patients' nutritional standard. The turning display allows easy reading from patient and user side. Scale operates by mains and batteries. Setting up of the scale is fast and straight forward. Supplied in just two parts and assembled by one special screw. An attachable telescopic height measuring rod is available as option.

Technical specification:

- Functions: Tare, hold, automatic zero setting, turning display, low battery indication, auto switch off, BMI calculation, kg/lbs, mother & child function
- Capacity & graduation: 250kg by 50g / 100g
- Dimensions & weight: 450 x 350 x 920mm 7,7kg
- Optionally: MZ10027-1 telescopic height measure

