

Column and floor scales



Practical and easy to use electronic floor scale for routine weight checks of patients, athletes or for any person concerned about their weight. The good-sized weighing surface will take small or large people giving them a very secure standing due to the tilde proof scale body. Whereas the scale has a full metal construction covered in a white ABS housing. The clear red dot matrix display allows the reading under all lighting conditions with a minimum of power consumption. The step on start together with an auto-hold function will give quick and reliable readings with a fine graduation of just 50 g. A non-slip robust rubber surface will add to the safe and secure sense when standing on the scale. At a weight of just under 1.6 kg the scale can easily be carried.

Technical specification:

Capacity:	200 kg /440 lbs
Graduation:	50 g / 0,1 lbs
Dimensions:	338 x 36 x 340 mm
Weight:	1,6 kg

Start via 'Step-on' or 'Tap-on' function, auto-hold function, automatic zero setting, auto switch off, slip-proof surface, kg/lbs switch-over. Power supply 4 x AAA batteries life time approx. 4000 cycles.